



Recommended Maximum Pacing Duration for Defibrillator Pads (In Hours)

40ms Pulse Width Example: **Philips** HeartStart MRx Monitor/Defibrillator/Pacemaker

		Pulse Amplitude (mA)																				
		10	20	30	40	50	60	70	80	90	100	110	120	130	140	150	160	170	180	190	200	
Pacing Rate (ppm)	30	24	24	24	22	20	16	14	12	10	9	8	8	7	6	5	4					
	40	24	24	24	16	15	14	12	10	9	8	7	6	5	5	4	4					
	50	24	24	22	14	14	12	10	9	8	7	6	6	5	4	4	3					
	60	24	24	20	12	12	10	10	9	8	6	6	5	5	4	3	3					
	70	24	24	16	10	10	9	8	8	6	5	4	4	3	3	3	3					
	80	24	24	14	8	8	7	6	6	4	4	3	3	2	2	2	2					
	90	24	24	12	8	7	6	5	5	4	3	3	2	2	2	2	1					
	100	24	24	10	6	6	5	5	4	4	3	3	2	2	2	2	1					
	110	24	24	8	6	5	4	4	3	3	3	2	2	2	2	1	1					
	120	24	24	8	7	4	4	3	3	3	2	2	2	2	1	1	1	1				
	130	24	24	7	6	4	4	3	3	2	2	2	2	1	1	1	1	1				
	140	24	24	6	5	3	3	3	2	2	2	2	1	1	1	1	1	1				
	150	24	24	6	4	3	3	2	2	2	2	2	1	1	1	1	1	1				
	160	24	24	5	4	3	2	2	2	2	2	2	1	1	1	1	1	1				
170	24	24	5	4	2	2	2	2	2	1	1	1	1	1	1	1	1					
180	24	24	4	3	2	2	2	2	2	1	1	1	1	1	1	1	1					