

Direct Connection for Your Specific Machine



Philips Medical



Physio-Control



Zoll Medical



Anderson

Defibrillator & Electrode Compatibility

Brand of Defibrillator	Model of Defibrillator	Medline Compatible Defibrillation Electrode Part Numbers
Philips Medical	» XLT » HeartStart XL » HP Codemaster XL (Plug Style Connector)	» MRX » MRX ALS » MDSPHL100 » MDSPHL200 » MDSPHL100LO » MDSPHL200LO » MDSPEPHL
	» HeartStart FR2+ » HeartStart Forerunner	» MDSPHL100 » MDSPHL200 » MDSPHL100LO » MDSPHL200LO
	» HeartStart FR3 » HeartStart FRx	Medline defibrillation electrodes are NOT compatible with these Philips devices
Physio-Control	» LP-9P » LP-10C » LP-12	» LP-15 » LP-20 » MDSPHY100 » MDSPHY200 » MDSPEDPHY
	» LP-500 » LP-1000	» LP-500 DPS » MDSPHY100 » MDSPHY200 » MDSPHY100LO » MDSPHY200LO
Zoll Medical	» M-Series » R-Series	» E-Series » X-Series » MDSZOL100 » MDSZOL200 » MDSPEZOL
	» Zoll AED PLUS » Zoll AED Pro	Medline defibrillation electrodes are NOT compatible with Zoll AED's.
Anderson	Compatible with defibrillators standardized to Anderson connector	» MDST100LOAC » MDST100AC » MDSPE100 » MDSC200LOAC » MDSC200AC

Medline Industries' pediatric defibrillation electrodes are NOT compatible with AED's. AED's cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

Contact your Medline Representative for any additional questions regarding electrode compatibility.

Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

FOLLOW US

<p>Medline United States 1-800-MEDLINE (633-5463) medline.com info@medline.com</p>	<p>Medline Canada 1-800-396-6996 medline.ca canada@medline.com</p>	<p>Medline México 01-800-831-0898 medlinemexico.com mexico@medline.com</p>
---	---	---

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.
©2018 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc.